

Identifying Age-Related Macular Degeneration from OCT Images

Jonathan Medley

Introduction

- Age-related macular degeneration (AMD):
 - Affects 8 million people in the United States
 - Early diagnosis leads to earlier treatment
 - Indicated by deposits called drusen
 - Symptoms – drusen structures
 - Diagnosis
 - Treatment


Introduction

- Stages of AMD:
 - 1: No AMD – a few small or no drusen
 - 2: Early AMD – several small, some intermediate
 - 3: Intermediate AMD – several intermediate 1+ large
 - 4: Advanced AMD – geographic atrophy and/or choroidal neovascularization – vision loss

Introduction

- Optical Coherence Tomography (OCT):
 - Uses interferometry to obtain images
 - Very high resolution of surface features

Introduction





<https://www.prescriber.co.uk/article/prevention-treatment-age-related-macular-degeneration/>

Introduction



- Dataset from Duke University
 - 384 patients
 - 269 AMD patients
 - 115 control
 - 100 images per patient
 - Depth varies with center of fovea at image 50

Original Images

AMD



Control




Methods




- Image filtering and binarization
- Isolating Retinal Pigment Epithelium
- Isolating RPEDC upper boundary as the working signal
- Filtering the signal
- Analyzing the power spectra and autocovariance of the signal

Processed Images

AMD



Control




Methods - Signal Acquisition




- Windowing / Thresholding – treat everything below 50% as black – “imadjust” function
- Fill Holes – Morphological operation
- Wiener Filter – adaptive low pass filter for smoothness – “wiener2” function
- Removing small image regions “bwareaopen” function
- Erosion and dilation – small disk – “imclose” function

Signals

AMD



Control




Methods – Signal Processing




- Removal of the mean
- Differentiation
- Median filtering
- Autocovariance calculation
- Power spectrum calculation
- Average power calculation
- Comparison of average Power for AMD vs Control Samples

Results - Autocovariances

AMD





Control





Results - Power Spectra


AMD



Control



Results - Comparisons



Summary

- RPEDC suggests useful clinical information for eye diseases
- OCT yields eye high quality images including RPEDC
- Processing OCT eye images provided a signal: RPEDC upper boundary
- Signal processing provided autocovariances and power spectra for analysis

Conclusions

- Autocovariance provided little information
- Average power was significantly higher in AMD samples than in control samples
- Greater number of samples need to be processed to create a classifier from the signal
- Average power may represent a useful method for classifying AMD from OCT images

References

- [1] Farsiu, S., Chiu, S. J., Oconnell, R. V., Folgar, F. A., Yuan, E., Izatt, J. A., & Toth, C. A. (2014). Quantitative Classification of Eyes with and without Intermediate Age-related Macular Degeneration Using Optical Coherence Tomography. *Ophthalmology*, 121(1), 162–172. doi: 10.1016/j.opthta.2013.07.013.