

# Identifying Age-Related Macular Degeneration from OCT Images

Jonathan Medley

# Background

- Age-related macular degeneration (AMD):
  - Affects 8 million people in the United States
  - Early diagnosis leads to earlier treatment
  - Indicated by deposits called drusen
  - Symptoms – drusen structures
  - Diagnosis
  - Treatment


# Background

- Stages of AMD:
  - 1: No AMD – a few small or no drusen
  - 2: Early AMD – several small, some intermediate
  - 3: Intermediate AMD – several intermediate 1+ large
  - 4: Advanced AMD – geographic atrophy and/or choroidal neovascularization – vision loss

# Background

- Optical Coherence Tomography (OCT):
  - Uses interferometry to obtain images
  - Very high resolution of surface features

# Background




<https://www.prescriber.co.uk/article/prevention-treatment-age-related-macular-degeneration/>





# Background

- Dataset from Duke University
  - 384 patients
  - 269 AMD patients
  - 115 control
  - 100 images per patient
  - Depth varies with center of fovea at image 50

# Sample Images – Normal







# Age-Related Macular Degeneration



# Methods

- Pre-filtering
- Binarization
- Filtering Binary Images
- Isolating Retinal Pigment Epithelium
- Comparing key measurements of the region
- Creating a threshold for classification

# Original Images




- Need filtering to remove noise, join segments, and isolate without losing information

# Pre-filtering





- Windowing – treat everything below 20% as black – “imadjust” function
- Wiener Filter – adaptive low pass filter for smoothness – “wiener2” function

# Pre-filtering



# Binarization


- Distinct regions already observed







# Filtering Binary Images

- Removing small objects – ignore regions smaller than 5000 pixels – “bwareaopen” function
- Erosion and dilation – small disk – “imopen” function
- Remove small objects again – 1000 pixels
- Fill holes in remaining objects – “imfill” function
- Dilation and erosion – slightly larger disk – “imclose” function
- Repeat

# Filtering Binary Images




# RPE Isolation




# RPE Parameters


- Convex Hull Area
- Perimeter
- Area
- Area/Perimeter
- Area/Major Axis




# Comparisons




# Comparisons




# Comparisons




# Comparisons



# Comparisons



# Comparisons



# Classification Threshold

- It was not possible to create a useable threshold from any of the variables analyzed in this model.
- Variations could be combined into a weighted system

# Conclusion

- Successful identification and isolation of RPE
- Consistent variations between variables of interest
- No single identifier to distinguish AMD and control patients
- Combination of variables may be possible for a classification system

# References

- [1] Farsiu, S., Chiu, S. J., O'Connell, R. V., Folgar, F. A., Yuan, E., Izatt, J. A., & Toth, C. A. (2014). Quantitative Classification of Eyes with and without Intermediate Age-related Macular Degeneration Using Optical Coherence Tomography. *Ophthalmology*, 121(1), 162–172. doi: 10.1016/j.ophtha.2013.07.013.